Anti peptidi Gliadina deamidati (DGP) IgG

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BIBLIOGRAFIA

1. Chartrand LJ, Seidman EG. Celiac disease is a lifelong disorder. Clin.Invest.Med., Vol. 19, 357-361, 1996
2. Cornell HJ. Coeliac disease: A review of the causative agents and their possible mechanisms of action. Amino Acids, Vol. 10, 1-19, 1996
3. Cronin CC, Feighery A, Ferriss JB, Liddy C, Shanahan F, Feighery C. High prevalence of celiac disease among patients with insulin-dependent (type I) diabetes mellitus. Am.J Gastroenterol., Vol. 92, 2210-2212, 1997
4. Jokinen J, Peters U, Maki M, Miettinen A, Collin P. Celiac sprue in patients with chronic oral mucosal symptoms. J Clin.Gastroenterol., Vol. 26, 23-26, 1998
5. Taminiau JA. Celiac disease. Curr.Opin.Pediatr., Vol. 8, 483-486, 1996
6. Williams CN. Celiac disease: past, present and future. Can.J Gastroenterol., Vol. 11, 647-649, 1997

Contatti:

InterMedical S.r.l. Via A.Genovesi,13 80010 Villaricca (NA) ITALY - Tel. +39 81 330 27 05 Fax +39 81 330 14 53

P. IVA 03426331215

e-mail *product specialist* :

mail@intermedical.it

INTERMEDICAL s.r.l.
Via A. Genovesi,13
80010 Villaricca(Na)-ITALY

Anti Deamidated Gliadin Peptide (DGP) IgG

Quantitative determination of IgG class antibodies against Deamidated Gliadin Peptide (DGP) in human serum or plasma

IVD

LOT

See external label

$\Sigma = 96$ tests

INTENDED USE

Anti Deamidated Gliadin Peptide (DGP) IgG is an indirect solid phase enzyme immunometric assay (ELISA) kit designed for the quantitative measurement of IgG class antibodies directed against Deamidated Gliadin Peptides (DGP) in human serum or plasma.

Anti-Deamidated Gliadin Peptide (DGP) IgG kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Celiac disease, also known as gluten sensitive enteropathy is primarily a disease of the infant organism. It is caused by a hypersensitivity reaction in response to gliadin, a protein being present in many cereals. This, non IgE mediated food allergy leads to massive malabsorption disturbances and is characterized by a complete atrophy of the villi and a hyperplasia of the crypts of the upper intestine.

Accordingly patients suffering from celiac disease must maintain a gluten free diet for the rest of their life.

Gliadins are proteins containing high amounts of the amino acids prolin and glutamine. These proteins belong to the nutritive tissue of the grain seeds of wheat, oat, barley and rye and are responsible for the baking properties of the flour.

Due to the possibilities of the highly specific and sensitive serological determination of IgA and IgG antibodies against DGP the invasive procedures of biopsies can be given up. In the past several biopsies have been done with patients when celiac disease was suspected, after a period of a gluten-free diet, and also after a specific gluten challenge; DGP antibodies titer has been proven to correlate very well with the morphological appearance of the mucosa of the upper intestine. It has been well documented that DGP antibodies level fall very quickly after a gluten free diet has begun and rise immediately after restoring gluten to the diet. Thus the serological test represents a reliable method to monitor patients, and in particular children and teenagers, for their adherence to the gluten-free diet.

2. PRINCIPLE

Anti-Deamidated Gliadin Peptide (DGP) IgG test is based on the binding of present antibodies in calibrators, controls or prediluted patient samples on the syntetic Deamidated Gliadin Peptides (DGP) coated on the inner surface of the wells. After a 30 minutes incubation the microplate is washed with wash buffer for removing non-reactive serum components.

An anti-human-IgG horseradish peroxidase conjugate solution recognizes IgG class antibodies bound to the immobilized antigens. After a 30 minutes incubation any excess enzyme conjugate, which is not specifically bound is washed away with wash buffer.

A chromogenic substrate solution containing TMB is dispensed into the wells. After 15 minutes of incubation the color development is stopped by adding the stop solution. The solutions color changes into yellow. The amount of color is directly proportional to the concentration of IgG antibodies present in the original sample.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

1. Anti-DGP Calibrators (5 vials, 1.2 mL each)

Phosphate buffer 0.1M, $\text{NaN}_3 < 0.1\%$, human serum

CAL0

CAL1

CAL2

CAL3

CAL4

2. Controls (2 vials, 1.2 mL each, ready to use)

Phosphate buffer 0.1M, $\text{NaN}_3 < 0.1\%$, human serum

Negative Control

Positive Control

3. Sample Diluent (1 vial, 100 mL)

Phosphate buffer 0.1 M, $\text{NaN}_3 < 0.1\%$

4. Conjugate (1 vial, 15 mL)

Anti h-IgG conjugated with horseradish peroxidise (HRP), BSA 0.1%, Proclin $< 0.0015\%$

5. DGP Coated Microplate

(1 breakable microplate coated with DGP)

6. TMB-Substrate (1 vial, 15 mL)

3,3',5,5'-tetramethylbenzidine 0.26 g/L, hydrogen peroxide 0.05%, Proclin $< 0.0015\%$

7. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15M

8. 10X Conc. Wash Solution (1 vial, 50 mL)

Phosphate buffer 0.2M, Proclin $< 0,0015\%$

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm)

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators and the Controls should be handled in the same manner as potentially infectious material.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of Sodium Azide (NaN₃) or Proclin 300^R as preservatives. Avoid the contact with skin or mucosa.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, allow scroll through large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use

disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinization and decontamination are effective in avoiding contamination of the conjugate; **this procedure is highly recommended when the kit is processed using analyzers which are not equipped with disposable tips.**

For this purpose, Intermedical supplies a separate decontamination reagent for cleaning needles.

- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibrators (C₀...C₄)

Since no international reference preparation for Anti-DGP antibodies is available, the assay system is calibrated in relative arbitrary units. The Calibrators are ready to use and have the following concentration:

	C ₀	C ₁	C ₂	C ₃	C ₄
AU/mL	0	15	30	60	240

Once opened, the Calibrators are stable 6 months at 2-8°C.

6.2. Preparation of the Sample

For determination of Anti-DGP human serum or plasma are the preferred sample matrixes.

All serum and plasma samples have to be prediluted with sample diluent 1:100; for example 10 µL of sample may be diluted with 990 µL of sample diluent.

The patients need not to be fasting, and no special preparations are necessary. Collect blood by venipuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation.

Samples may be stored refrigerated at 2-8°C for at least 5 days. For longer storage of up to six months samples should be stored frozen at -20°C. To avoid repeated thawing and freezing the samples should be aliquoted.

Neither Bilirubin nor Hemolysis have significant effect on the procedure.

The Controls are ready to use.

6.3. Preparation of the Wash Solution

Dilute the contents of each vial of the buffered wash solution concentrate (10X) with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C. In concentrated wash solution is possible to observe the presence of crystals; in this case mix at room temperature until the complete dissolution of crystals; for greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care to transfer completely the crystals, then mix until crystals are completely dissolved.

6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.**
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₄), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Controls	Blank
Calibrator C ₀ -C ₄	100 µL		
Controls		100 µL	
Diluted Sample		100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 µL of diluted wash solution			
Conjugate	100 µL	100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 µL of diluted wash solution.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate 15 minutes in the dark at room temperature (22-28°C).			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against Blank within 5 minutes.			

7. RESULTS

7.1. Calibration curve

For Anti-DGP IgG a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice. Smoothed-Spline Approximation and log-log coordinates are also suitable. However we recommend using a Lin-Log curve.

First calculate the averaged optical densities for each calibrator well. Use lin-log graph paper and plot the averaged optical density of each calibrator versus the concentration. Draw the best fitting curve approximating the path of all calibrator points. The calibrator points may also be connected with straight line segments. The concentration of unknowns may then be estimated from the calibration curve by interpolation.

Typical Results (example only)

The figures below show typical results for Anti- DGP IgG. These data are intended for illustration only and should not be used to calculate results from another run.

N	OD1	OD2	mean	C1	C2	mean	CV%
CAL0	0,011	0,011	0,011	0,09	0,09	0,09	6E-7
CAL1	0,165	0,161	0,163	14,74	14,36	14,55	1,89
CAL2	0,322	0,324	0,323	30,46	30,67	30,57	0,48
CAL3	0,590	0,590	0,590	59,78	59,78	59,78	4E-7
CAL4	1,694	1,768	1,731	232,1	248,1	240,1	4,71

8. REFERENCE VALUES

In a normal range study with serum samples from healthy blood donors the following ranges have been established with the Anti-DGP tests:

	Anti DGP IgG (AU/mL)
Negative	< 15
Equivocal	15 - 30
Positive	> 30

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken individually.

It is recommended that each laboratory establishes its own normal and pathological ranges of serum Anti-DGP.

8.1. Specificity

Comparison test against a commercial reference kit, performed on 63 sera (32 of them positive sera and 31 negative sera) showed a 100.0% specificity.

8.2. Sensitivity

Comparison test against a commercial reference kit, performed on 63 sera (32 of them positive sera and 31 negative sera) showed a 94.1% sensitivity.

8.3. Detection limit

The lowest concentration of anti-DGP IgG that can be distinguished from the Calibrator 0 is 0.13 AU/mL with a confidence limit of 95%.

8.4. Precision and reproducibility

8.4.1. Intra-Assay

Within run variation was determined by replicate 16 times two different sera with values in the range of calibration curve. The within assay variability is \leq 3.8%

8.4.2. Inter-Assay

Between run variation was determined by replicate the measurements of two different control sera with different lots of kits and/or different mix of lots of reagents. The between assay variability is \leq 7.8%.

9. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

1. Chartrand LJ, Seidman EG. Celiac disease is a lifelong disorder. Clin.Invest.Med., Vol. 19, 357-361, 1996
2. Cornell HJ. Coeliac disease: A review of the causative agents and their possible mechanisms of action. Amino Acids, Vol. 10, 1-19, 1996
3. Cronin CC, Feighery A, Ferriss JB, Liddy C, Shanahan F, Feighery C. High prevalence of celiac disease among patients with insulin-dependent (type I) diabetes mellitus. Am.J Gastroenterol., Vol. 92, 2210-2212, 1997
4. Jokinen J, Peters U, Maki M, Miettinen A, Collin P. Celiac sprue in patients with chronic oral mucosal symptoms. J Clin.Gastroenterol., Vol. 26, 23-26, 1998
5. Taminiu JA. Celiac disease. Curr.Opin.Pediatr., Vol. 8, 483-486, 1996
6. Williams CN. Celiac disease: past, present and future. Can.J Gastroenterol., Vol. 11, 647-649, 1997

Contact:

InterMedical S.r.l. Via A.Genovesi,13 80010 Villaricca (NA) ITALY - Tel. +39 81 330 27 05 Fax +39 81 330 14 53

P. IVA 03426331215

e-mail *product specialist* :

mail@intermedical.it

INTERMEDICAL s.r.l.
Via A. Genovesi,13
80010 Villaricca(Na)-ITALY

	Intermedical	
	PACKAGING INFORMATION SHEET	

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


GB
Explanation of symbols

FR
Explication des symboles

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DE
Verwendete Symbole

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	DE In vitro Diagnostikum ES Producto sanitario para diagnóstico In vitro FR Dispositif medical de diagnostic in vitro GB In vitro Diagnostic Medical Device IT Dispositivo medico-diagnostico in vitro PT Dispositivos medicos de diagnostico in vitro		DE Hergestellt von ES Elaborado por FR Fabriqué par GB Manufacturer IT Produttore PT Produzido por
REF	DE Bestellnummer ES Número de catálogo FR Références du catalogue GB Catalogue number IT Numero di Catalogo PT Número do catálogo	 yyyy-mm	DE Herstellungs datum ES Fecha de fabricacion FR Date de fabrication GB Date of manufacture IT Data di produzione PT Data de produção
 yyyy-mm-dd	DE Verwendbar bis ES Estable hasta (usar antes de último día del mes) FR Utiliser avant (dernier jour du mois indiqué) GB Use by (last day of the month) IT Utilizzare prima del (ultimo giorno del mese) PT Utilizar (antes ultimo dia do mês)		DE Biogefährdung ES Riesco biológico FR Risque biologique GB Biological risk IT Rischio biologico PT Risco biológico
 	DE Gebrauchsanweisung beachten ES Consultar las instrucciones FR Consulter le mode d'emploi GB Consult instructions for use IT Consultare le istruzioni per l'uso PT Consultar instruções para uso		DE Chargenbezeichnung ES Codigo de lote FR Numero de lot GB Batch code IT Codice del lotto PT Codigo do lote
 Σ = xx	DE Ausreichend für "n" Tests ES Contenido suficiente para "n" tests FR Contenu suffisant pour "n" tests GB Contains sufficient for "n" tests IT Contenuto sufficiente per "n" saggi PT Contém o suficiente para "n" testes		DE Inhalt ES Contenido del estuche FR Contenu du coffret GB Contents of kit IT Contenuto del kit PT Conteúdo do kit
 Max Min	DE Temperaturbereich ES Límitación de temperatura FR Limites de température de conservation GB Temperature limitation IT Limiti di temperatura PT Temperaturas limites de conservação		

	Intermedical	
	PACKAGING INFORMATION SHEET	

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ERROR POSSIBLE CAUSES / SUGGESTIONS

No colorimetric reaction

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation